

Identifying Ventilator Dyssynchronies via Waveform Assessment

Tim France BS, RRT

Hamilton Medical Inc.

Clinical Applications Specialist, USA

Disclosures

- Currently an employee of Hamilton Medical Inc.

Objectives

- Explain Importance of ventilator synchrony
- Identify common synchrony issues via waveform
- Define Synchrony vs Dyssynchrony
- Learn how to make basic vent changes to marry the vent to the patient

Why is synchrony important

- Patient comfort
- Less sedation
- Less time on the vent
- Less alarms

Definitions for ventilator/patient interactions

- Synchrony- vent and patient are in phase working together
- Dyssynchrony- vent and patient not in phase. Typically vent is not responding to changes in patient conditions

Pulmonary Mechanics

- P_{0.1}- Pressure at the first .1 seconds of a breath. Measures a patient's central drive to breath. The more negative the number the higher the drive.
- Compliance
 - Chest wall vs Lung
- Resistance
 - Inspiratory- measures force prior to lower airways (Asthma)
 - Expiratory- measure force during exhalation. Collapsed airways (COPD)
- RC_{exp} Expiratory time constants
 - Longer for COPD
 - Shorter for ARDS



Missed inspiratory attempts

Assess

- P0.1
- Pt effort
- Trigger setting
- Autopeep

Intervention

- ↓ trigger setting
- ↑ PEEP to match or decrease VE



Inadvertent Triggering (leak)

Assess

- Patient effort, check P0.1, circuit integrity

Intervention

- Increase Trigger setting
- Fix leak



Inadvertent Triggering

- Assessment
 - Assess drive, check P0.1,
- Causes
 - Heart Rate
 - External Device
 - Pacemaker
- Intervention
 - Increase trigger setting
 - Manipulate external device



CMV No patient effort



Flow Starvation in Flow Controlled modes

- Assess P0.1, Pminimum
- Actions
 - ↑ PF
 - Switch to Flow variable mode (PCV, Spont)



Increased Patient work in Pressure Targeted modes

- Assess P0.1
- Actions
 - ↑ support until flow waveform decelerates if significant patient work is observed



Increased Patient work Spontaneous mode

- Check flow waveform for full deceleration
- Assess P0.1
- Actions
 - ↑ support until flow waveform decelerates if significant patient work is observed



CMV No patient effort



P_{aw} Increase because of compliance changes

Paw increase at end inspiration

- ARDS
- Fibrotic
- Chest wall compliance issues

Actions

- Increase PEEP (if early ARDS or chest wall issues)
- Check Tidal volume




Plateau and Peak are the same

- Perform inspiratory hold



P_{aw} increase because of airway resistance

- Paw increase is in the beginning of Inspiration
- Causes:
 - COPD, Asthma, small ET tube
- Interventions
 - Bronchodilator
 - Heliox
 - Consider larger tube
 - Increase PEEP



Plateau with increase in Resistance

- Plateau is significantly less than Paw



ETS (Expiratory Trigger Sensitivity)

- ETS 5%
- Once flow decelerates to 5% of initial flow exhalation starts
- If PF is 100 lpm once flow decelerates to 5 lpm exhalation starts
- I-time is longer
- Vt is larger



ETS 70%

- Once flow decelerates to 70% of initial flow exhalation starts.
- I-time is shorter
- Vt is lower



Pramp (Pressure Ramp)

- How fast pressure setting is reached
- Pramp of 5 ms
 - Initial flow is higher
 - Vt is higher
- Pramp of 200
 - Initial flow is lower
 - Vt is lower
 - Paw max is delayed



Pramp to Aggressive

Assess

- Look for spike at beginning of inspiration
- Patient comfort
- ET tube size

Intervention

- Increase Pramp setting till spike goes away



Late cycling

- Assess

- I-time setting
- ETS if in spont mode

- Intervention

- Decrease I-time till bump goes away
- Increase ETS setting
- Place in spont mode to assess spontaneous I-time



Early cycling

- **Assess**

- I-time
- ETS if in spontaneous mode

- **Intervention**

- Increase I-time till double breaths stop
- Decrease ETS %
- Place patient in spont mode to find spontaneous I-time



Oscillations in the respiratory circuitry

Causes

- Secretions
- Condensate

Intervention

- Suction
- Clear Circuit

Conclusion

- Assess waveforms for synchrony
- Assess waveforms to tweak vent settings

Thank You For Your Time

Tim France BS, RRT

Clinical Applications Specialist

Tim.France@hamiltonmedical.net

Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

☎ +41 58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com